



AirForce Mini

Compact Compressor Nebulizer

TECHNICAL DATA

Dimensions

Length: 111 mm

Height: 45 mm

Width: 68 mm

Weight

217 grams

Particle Size:

MMAD 4.12 micrometres

Nebulization rate: 0.25ml/min.

Min. filling volume: 2 ml

Max. filling volume: 6 ml

Residual volume: 0.8 ml

The AirForce Mini from DriveDeVilbiss Healthcare is a battery-operated portable inhalation therapy system. Whether in the office, at school or on the road, the AirForce Mini gives the user complete independence.

SPECIFICATIONS

- Lightweight
- Low noise
- Integrated battery with LED-status display
- Stylish design
- User friendly and easy to clean
- One of the smallest nebulizer systems in the market

Designed using modern day compressor-technology, the AirForce Mini produces a high percentage of respirable aerosol particles suitable for treating respiratory conditions of the upper and lower airways. Exceptional performance and value; the AirForce Mini is the ideal compact solution for aerosol therapy which is suitable for both adults and children.

Contents included:

AirForce Mini compressor nebulizer, mouthpiece, tubing, replacement filters (5x), adults mask, pediatric mask, nosepiece, Micro-USB-cable, nebulizer, instruction guide, bag and power adapter.

Drive DeVilbiss International

DeVilbiss Healthcare GmbH
Kamenzer Straße 3
68309 Mannheim, Germany

UK +44 (0) 121 / 521 31 40
DE +49 (0) 621 / 1 78 98 0
kontakt@devilbisshc.com

www.drivedevilbiss-int.com

drive
DeVilbiss[®]
HEALTHCARE