

VacuAide



DOEVLBISS

Suction



Portable Suction Unit

VacuAide 7305

The 7305 unit is ideal for home use and meets the stringent ISO standards for emergency transport suction. The 7305 is small and lightweight, but with a high performance capability that meets the high demand needs of the patient. The 7305 is available in two models, the high performance 7305 D-U or the 7305 P-U with added rechargeable battery and an optional 12 V vehicle charger making the 7305 P-U a truly portable unit.

Features and benefits (both versions)

- Switch mode power supply allowing operation on any AC voltage 100 V - 240 V
- High-performance capabilities, up to 550 mmHg vacuum and 27 l/min flow
- Small, lightweight design for ease of portability
- Integrated bacterial filter
- Automatic float shut-off
- 800 ml disposable collection bottle
- Optional 1200 ml re-usable collection bottle

Battery version 7305 P-U

- Rechargeable battery with 60 min. run time
- Automatically recharges when connected to the power supply
- Canvas Carrying Case with shoulder strap included
- 12 V DC capability for vehicle charging

Areas of Application

- Bronchial / Oral Hygiene
- A means to provide effective airways clearance
- A means to control infection



SUNRISE
MEDICAL

SPECIFICATIONS

Dimensions (HxWxD)	23 x 18 x 20 cm	AC Operation	110 - 240 V; 50/60 Hz
Weight 7305 D-U	1.7 kg (3.75 lb) ¹⁾	Internal Battery runtime	60 min. (7305 P-U)
Weight 7305 P-U	2.9 kg (6.40 lb) ¹⁾	¹⁾ Weight of the unit does not include the AC adapter.	
Vacuum Range	80 - 550 mmHg ²⁾	²⁾ Conditions may vary based on altitude above sea level, barometric pressure changes and system configuration	
Flow Range	0.0 - 27 l/min ²⁾		

Sunrise Medical has a policy of continuing improvement, we therefore reserve the right to vary details at any time without notice.

Distributor:

0804/2000/ST-serv/0690137_UK

Sunrise Medical Limited
High Street Wollaston
West Midlands DY8 4PS
England
Phone + 44 (0) 13 84 44 68 58
Fax + 44 (0) 13 84 44 66 28
www.sunrisemedical.co.uk
devilbiss@sunmed.co.uk

CE 0044

 **SUNRISE
MEDICAL**